

DMS Series: Goniometric Display Testing



- » Angular and areal 2D dependency of luminance, color and contrast
- » 6 motorized axes for angular, directional and lateral scanning
- » Electro-optical transfer characteristics (EOTF) and flicker
- » Spectral diffuse/direct reflection and transmission
- » 300 μm sphere of confusion (6-axis movement)
- » DUT size up to 72" (~ 2m)
- » Positioning resolution 0.01° / 0.01 mm
- » Thermal chamber for -40 °C to +105 °C



**Large global customer base,
worldwide service!**