

LumiCam 2400B Series

2D Imaging Colorimeter

Key features at a glance

- ▲ Motorized objective lenses
- ▲ Up to 30 % shorter measurement time due to optimized stepper motor actuation
- ▲ Approx. 20 % smaller footprint for compact design
- ▲ Various analysis features such as spotmeter, polygons, and flexible polylines
- ▲ Straightforward determination of luminance distributions, color coordinates / temperature, Planck distance, dominant wavelength, or contrast



\\ TECHNICAL SPECIFICATIONS

LumiCam 2400B	Mono	Color	Advanced
Quantities			
Photometric quantities	Luminance (cd/m ²), luminous intensity (cd), contrast		
Colorimetric quantities	-	Color coordinates (x,y), color coordinates (u',v'), tristimulus values (X, Y, Z), dominant wavelength (nm), color saturation, correlated color temperature CCT (K)	
Camera data			
Sensor	Sony IMX264LLR CMOS Sensor		
Sensor size	2/3", 11.1 mm diagonal		
Effective number of pixels (h x v)	2428 x 2028 (5 MP)		
Pixel size	3.45 μm x 3.45 μm		
AD converter	12 bit		
Exposure time	40 μs to 30 s		
Luminance measurement			
Measurement range ¹⁾	0.3 mcd/m ² – 2.5 Mcd/m ²		
Extended measurement range ²⁾	2.5 x 10 ¹⁰ cd/m ²		
Measurement time incl. data saving time (at 10 cd/m ²) ³⁾	1.9 s		
Measurement time incl. data saving time (at 100 cd/m ²) ³⁾	1.8 s		
Accuracy for std. illuminant A ⁴⁾	±3 %	±3 %	±3 %
Accuracy for LED color light ⁸⁾	-	-	±2 %
Repeatability ⁵⁾	±0.03 %		
Linearity	±0.5 %		
Uniformity ⁶⁾	±0.5 %		
Filter match ⁷⁾	f ₁ ¹ < 3 %		

\\ TECHNICAL SPECIFICATIONS

Color measurement			
Measurement time incl. data saving time (at 10 cd/m ²) ³⁾	-	10 s	16 s
Measurement time incl. data saving time (at 100 cd/m ²) ³⁾	-	8 s	12 s
Accuracy (x, y) for std. illuminant A ⁴⁾	-	±0.003	±0.003
Accuracy (x, y) for color light ⁸⁾	-	±0.010	±0.010
Accuracy (x, y) for LED color light ⁹⁾	-	-	±0.005
Repeatability (x, y) ⁵⁾	-	±0.0001	±0.0003
Uniformity (x, y) ⁶⁾	-	±0.001	±0.001
Filter match	f _i ' (Y) < 3 %	f _i ' (Xb) < 6 % f _i ' (Xr) < 6 % f _i ' (Y) < 3 % f _i ' (Z) < 4 %	f _i ' (Xb) < 6 % f _i ' (Xr) < 6 % f _i ' (Y) < 3 % f _i ' (Z) < 4 %
General			
Interface	Gigabit Ethernet		
Operating system	Windows 7 (32/64 bit), Windows 10 (64 bit)		
Dimensions (L x W x H) (incl. 50 mm objective lens, no handle)	210 mm x 105 mm x 100 mm	235 mm x 119 mm x 133 mm	
Weight	1.5 kg	3 kg	
Power supply	24 V external		
Operating conditions	10 to 40 °C, max. 70 % relative humidity (non-condensing)		

Instrument Systems is continually working to develop and improve products. Technical changes, errors or misprints do not constitute grounds for compensation. The company's terms of delivery and payment apply in all other respects.

¹⁾ The lower limit corresponds to the maximal exposure time and smallest f-number at SNR 10:1 and vice versa.
²⁾ Valid for OD4 filter.
³⁾ Value calculated from 100 repetitions taking data saving time into account.
⁴⁾ Calculated from 100 repetitions; refers to the deviation of the mean from the reference value.
⁵⁾ Calculated from 50 repetitions. Refers to the double standard deviation of the measured values.

⁶⁾ Refers to the maximum deviation from the mean value calculated from flat-field verification image.
⁷⁾ Deviation of the filter transmission from the V(λ) curve integrated over the entire visible spectrum.
⁸⁾ Maximum deviation from the reference source (illuminant A with set of color glass filters).
⁹⁾ Derived from 20 repetitions for R, G, and B - LEDs with optimized signal level. Refers to the double standard deviation.