

CAS 140D

Change Notice:

CAS 140D Series Design Update

Date of notice	Valid from
2026-07-15	All orders placed after 15th September 2026
Description of change	
For all purchase orders placed after 15th September 2026, the CAS 140D will feature mechanical design modifications to the front switch and the fixation of the interface module:	
1. Omission of the front switch	
New: without front switch 	Old: with front switch 
2. Replacement of captive screws with standard screws for the interface module	
New: with standard screws 	Old: with captive screws 
Reason for change	
The design adjustments reflect common usage behavior and ensure stable long-term availability of all components without affecting the device's optical or electrical performance.	

Impact on device performance

The mechanical design changes have no impact on the optical or electrical performance of the CAS 140D. All specifications such as spectral capture, integration times and straylight performance remain fully unchanged.

The modifications are limited to two mechanical aspects:

1. Omission of the front switch

The spectrometer is now switched off exclusively via the rear power switch, ensuring consistent operation across all configurations.

2. Use of standard screws for the interface module

The interface module can be replaced using standard screws. When doing so, please ensure that the screws used for fixation are retained.

Affected Products/part numbers:

CAS140D-23594	CAS140D151U2L	CAS140D154E1B	CAS140D156U2D
CAS140D0941U2K1	CAS140D151U3A	CAS140D154E2B	CAS140D156U2F
CAS140D1130U1Q3	CAS140D151U3C	CAS140D154P1H	CAS140D156U2H
CAS140D151E1A	CAS140D151U3E	CAS140D154U1B	CAS140D156U3B
CAS140D151E1C	CAS140D151U3L	CAS140D154U1F	CAS140D157E1B
CAS140D151E1E	CAS140D152U1D	CAS140D154U2B	CAS140D157E1F
CAS140D151E1L	CAS140D153E1A	CAS140D154U2F	CAS140D157E1H
CAS140D151E2E	CAS140D153E1C	CAS140D154U3B	CAS140D157E3F
CAS140D151E2L	CAS140D153E2A	CAS140D156E1B	CAS140D157P1B
CAS140D151E3C	CAS140D153E3A	CAS140D156E2B	CAS140D157P1H
CAS140D151E3E	CAS140D153P1A	CAS140D156E2F	CAS140D157U1B
CAS140D151E3L	CAS140D153P1C	CAS140D156E3F	CAS140D157U1D
CAS140D151P1A	CAS140D153P2C	CAS140D156P1B	CAS140D157U1F
CAS140D151P1L	CAS140D153P2G	CAS140D156P1D	CAS140D157U2B
CAS140D151P2L	CAS140D153P3A	CAS140D156P1F	CAS140D157U2F
CAS140D151P3L	CAS140D153P3E	CAS140D156P1H	CAS140D157U3B
CAS140D151U1A	CAS140D153U1A	CAS140D156P2D	CAS140D171U1A
CAS140D151U1C	CAS140D153U1E	CAS140D156P2H	CAS140D171U1F
CAS140D151U1E	CAS140D153U1I	CAS140D156P3F	CAS140D171U2A
CAS140D151U1L	CAS140D153U2A	CAS140D156U1B	
CAS140D151U1O	CAS140D153U3A	CAS140D156U1D	
CAS140D151U2A	CAS140D153U3C	CAS140D156U1F	
CAS140D151U2C	CAS140D153U3I	CAS140D156U2B	

— All SETs containing the above part numbers —